

KX5600

Full digital B mode ultrasound scanner

NEW



Net weight of main unit is 4.5Kg, smart and light, convenient for carrying.



Large-capacity lithium battery can drive the main unit to work for more than 200 minutes.




The hard disk can be extended to 4G to get 5000 frames permanent image storage.

KAI XIN

XUZHOU KAI XIN ELECTRONIC INSTRUMENT CO., LTD.

Humanized design

- Rapidly one-key operation includes one-key storage, one-key browse and one-key print.
- Unique CPU temperature control function, when temperature is too high, it can automatically shut down.
- Rotatable LED Panel to make the users get the best image visual sense from different angles.
- One-key to light the keyboard, suitable for operation in dark room.
- Practical function knob design, use A key to realize quick adjustment in a variety of conditions.




Powerful data management system

- Optimization of the query image function, can realize 1-2 seconds response to your needs.
- Neoteric internal workstation multi-plan template, quickly finish data, report manage function, realize process multi medical histories function.

Abundant software package

- Equipped with professional gravel package which can realize real-time positioning measurement.
- OB measurement and analysis: 15 kinds OB table and several GA measurement, AFI, Fetus growth curve and OB report.
- GYN measurement and analysis: general GYN software packages and GYN report.
- Urological measurement and analysis: residual urine, volume of prostate, Urological report.



High quality clinical image

- High resolution and full screen display, meet your vision needs.
- Adjustable LED brightness, convenient and practical.



Standard configuration:

- Main unit: 1 pc
- 3.5MHz multi-frequency convex probe: 1 pc

Optional configuration:

- 6.5MHz multi-frequency trans-vaginal probe
- 7.5MHz multi-frequency linear probe
- Lithium Battery
- Charger
- Trolley

ADD: Kaixin Mansion, C-01, Economic Development Zone, Xuzhou, Jiangsu, China.

Zip Code: 221004

Tel: +86-516-87732932 / 87733758

Fax: +86-516-87732932 / 87792848

Http://www.kxele.com

E-mail: info@kxele.com / service@kxele.com

KAI XIN
XUZHOU KAI XIN ELECTRONIC INSTRUMENT CO., LTD.

THE PRODUCT TAKES REAL OBJECT AS QUASIF ANY FLUCTUATION HAPPENS, FORGIVING NO ANOTHER NOTICE IN THE LINE.